

### **TECHNICAL SPECIFICATIONS**

### **GENERAL CONSTRUCTION FEATURES**

- battery and mains operation
- stable base construction
- PU material column casing, easily removable
- base cover made of ABS plastic or stainless steel, both resistant to impact, breakage and disinfectant agents, easy to clean
- mobile with four double-swivel antistatic castors plus fifth wheel (ø100 mm) adjusted and locked electrically via hand switch
- · manual SOS castors release
- tabletop is divided into 4/5 sections
- · fully electric operation via hand switch
- · headrest adapter with 3 attachment options
- stainless steel fixation rails for accessories attachment (25x10mm)

#### **ELECTRICAL DATA**


- · features two batteries
- batteries recharged from the mains power supply; 100-240 V, 50-60 Hz, I: Max 4,0 A
- battery charge level with optical display on the hand switch
- Safety class I type B; the enclosure leakage current meets the requirements of the patient leakage current for CF conditions as per EN 60601-1




#### **SMART CONTROL FEATURES**

- memory function with 2 memory positions
- 2 pre-set positions: Trendeleburg, zero
- optical display for battery charge level
- remote access via Bluetooth for service diagnostics

# **OPERATING TABLE Y, ORTY 2400**



## **OPERATING TABLE Y, ORTY 2500**



TECHNICAL DATA	ORTY-2400	ORTY-2500
tabletop	4-section	5-section
length of tabletop incl. headrest	1940 - 2040 mm	1940 - 2040 mm
length of tabletop without headrest	1610 mm	1610 mm
narrow headrest (W x L)	300 x 330 mm	300 x 330 mm
wide headrest (W x L)	550 x 330 mm	550 x 330 mm
removable sections	head	head & leg
length of back section	550 mm	550 mm
length of middle section	540 mm	540 mm
length of seat section	/	260 mm
length of leg section	520 mm	/
length of detachable leg section	/	260 mm
width of tabletop	550 mm	550 mm
width with fixation rails	600 mm	600 mm
safe working load	320 kg	320 kg
maximum patient weight	230 kg	230 kg
total weight	145 kg	147 kg
thickness of upholstery pads	80/100 mm	80/100 mm



ADJUSTMENT RANGE	ORTY-2400	ORTY-2500
height adjustment range (without padding)	580 - 980 mm	580 - 980 mm
Trendelenburg/Anti-Trendelenburg	± 20°	± 20°
lateral tilt	± 20°	± 20°
back section	+73°	+73°
leg section	-90°	1
seat section with detachable leg section	/	-90°

